


GVP



1300 838 438 – www.generalvetproducts.com.au

NUWEIGH'S Veterinary Scales with large stainless-steel platform uses four precision loadcells for high accuracy readings.



4x
LOAD
CELLS




FUNCTIONS

- HOLD
- TARE
- UNIT
- LCD

FEATURES

- Four precision loadcells
- Auto OFF function for power saving
- Animal hold function for effortless recordings
- Adjustable feet
- Reads in both kgs, oz and lbs
- 9V power adaptor or 4 x AA battery

DIMENSIONS



Capacity	Divisions	Part
250kg	100g	No.

#36617

1300 838 438 / 1300 VET GET

sales@generalvetproducts.com.au

www.generalvetproducts.com.au