

# iM3 - The Veterinary Dental Company

## iM3 Pro 2000 LED Dental Unit

The iM3 Pro 2000 combines all the features iM3 have to offer in a compact self-contained dental system. The added storage under the Corian® top and the option of having two LED high speed handpieces make this system a must for the busy veterinary practice.



Magnetic bur holder to prevent burs falling on the ground.

Corian® top, impervious to bacteria, easy to clean and maintain with storage underneath

Autoclavable, adjustable speed suction

The Pro 2000 is supplied with an iM3 Advantage Swivel LED High Speed Handpiece. The 360-degree swivel allows greater control of the handpiece and reduces hand fatigue.

Shown with optional second high speed handpiece

Height adjustable, stainless steel stand

Large air tank to prevent the compressor motor overheating

Quiet, low maintenance oil-free compressor with CE approval

Large hospital grade castors

“Our technicians fought over who got to use the iM3 so we had to buy a second one!”

PRO2000 User, Oxford, OH, USA

**iM3 Pro 2000 LED with Oil-free Compressor**

Cat. No  
U3080

Don't forget to add iM3 Straw

Cat. No  
PS001



# iM3 - The Veterinary Dental Company



Shown with optional second high speed hand piece



## Features and Benefits of the iM3 Pro 2000 Dental Unit

The Pro 2000 is fitted with an iM3 Advantage Swivel LED High Speed Handpiece that features 360 degree swivel reducing unnecessary hand fatigue and allowing greater control of the handpiece.

The Pro 2000 comes with the following quality handpieces:

Variable suction – Three-way syringe (air, water and mist) – Advantage Low Speed Handpiece – Advantage Swivel LED High Speed Handpiece – LED light wand – 4:1 Advantage Swivel Straight Nose Cone for polishing – 1:1 Advantage Straight Nose Cone for polishing – 1:1 Straight Nose Cone for restorative dental attachments as well as rabbit and orthopaedic surgery using HP burs.

All iM3 low and high speed handpieces are autoclavable.

Furthermore, the Pro 2000 dental unit is supplied with an additional handpiece cable which gives you the option of having two LED high speed handpieces or an additional low speed handpiece with water cooling.

This high end dental unit is recommended for the busy veterinary practice and it gives the practice all the features iM3 have to offer in a compact self-contained dental system.

The Pro 2000 offers a large useable storage space. The smart hinge system means the Corian® top stays flat when its top is opened. Unlike other dental machines with tops that open up from a hinge at the back, the iM3 hinge system means the items on top do not fall on the floor.

The built-in magnetic bur holder prevents the burs from falling on the floor.

The Pro 2000 includes an LED light wand for oral illumination. The light allows you to explore the mouth and see cracks in teeth, while providing additional light in the mouth.

The Pro 2000 is fitted with a two bottle water system with auto pressure release. This allows the user to choose between CLS (chlorhexidine based solution) and distilled water depending on the procedure.

It's difficult and expensive to send a large dental machine away for servicing. The Pro 2000 is built to last. The cartridge and water valves designed by iM3 will keep working well beyond the life of an ordinary valve. Our systems are unique, and our engineering is so functional, that you will be impressed by the minimal amount of maintenance required to keep your system in top working order.

### **iM3 Pro 2000 LED with Oil-free Compressor**

Cat. No  
U3080

### **iM3 Pro 2000 LED SW**

(Slow speed with water cooling, not available in all countries)

U3080SW

Don't forget to add iM3 Straw

Cat. No  
PS001