



GVP Pty Ltd



Starting Our

**END OF YEAR
SALE**

Early



**OFFERS VALID UNTIL 30/06/2024
OR UNTIL STOCKS LAST
ORDER AND PAY NOW TO SECURE
THESE PRICES FOR THIS
FINANCIAL YEAR.**

**WE CAN HOLD YOUR GOODS UNTIL
YOUR READY FOR DELIVERY** 🚚

GVP
General Vet Products

EOFY Starting Early



OFFERS VALID UNTIL 30/06/2024 OR UNTIL STOCKS LAST
ORDER AND PAY NOW TO SECURE THESE PRICES FOR THIS FINANCIAL YEAR.

WE CAN HOLD YOUR GOODS UNTIL YOUR READY FOR DELIVERY 🚚



RRP \$990

(Excludes GST and Freight)

**Save even more
with our
EOFY Sale**

**Login in to see
how much cheaper
it is for GVP
Members**

Electric Low Lifting Exam Table - GVPTE-FT-801

Includes;

- Power Points
- Tools organising drawer
- Convenient Hanging hooks
- **Powered by Denmark Linak Actuator & Transformer**
- Maximum load capacity of 150kg (300lbs)
- Size: 125cm(h) x 65cm(w) x 25-95cm(adjustable height range)



GVP
General Vet Products

EOFY

SALE!!

Starting Early

AT GENERAL VET PRODUCTS

VetSpecs® Capno-5

Bluetooth Micro-flow™ Side-stream Capnography Monitor

GVPVSP-CAPNO5

Designed for monitoring the ventilation of intubated or non-intubated patients in surgery, dentistry, procedural sedation, pre-op preparation, anaesthetic recovery, emergency, critical care, and transport, the **VetSpecs® Capno-5 provides**

- Real-time CO2 Waveforms
- End-tidal CO2 Reading (EtCO2)
- Respiratory Rate (RR)
- Fractional Inspired CO2 (FiCO2)
- EtCO2 and RR Trends
- EtCO2 and Apnea Alarms



Includes:

- Vetspecs Handheld EtCO2 Monitor - SideStream
- Rechargeable Li-ion battery
- Protective cover
- Recharge normal power outlet or USB

BUY 2 & GET 20% OFF

(Excludes GST and Freight)

GVP
General Vet Products

EOFY

SALE!!

Starting Early

At General Vet Products

NORMAL PRICE ~~\$865.00~~

GVP MEMBER PRICE \$821.75

BONUS

**BUY 2 GET THE
3RD HALF PRICE!**

(Available until stocks last | Excludes GST and Freight)

The **800I IV Infusion Pump** is an easy to use, yet sophisticated pump that gives vet teams the freedom to use an extensive range of standard Giving Sets.

- Small, lightweight & compact
- Patented "Double Pulse" compensation technique ensures good linearity even at low delivery rates.
- Double CPU system makes the running status super safe. Delivery will be cut off, even if control circuit is powered off, and machine will alarm.
- Rechargeable Back-up Battery supports at least 3 hours continuous operation.
- Different fluid delivery modes available: Fluid delivery rate, Delivery time rate, Dripping Rate.
- Calibrate in-house with ease; no dedicated brand giving set needed



GVPPM-800I

GVP
General Vet Products

EOFY

SALE!



MIRUSTM

Linear Stapler

SAVE
OVER 20% OFF
PROMO Pack

MIRUSTM - The Symbol of Excellence which ensures superior surgical outcomes.

MIRUS Reloadable Linear Stapler - Gives two staggered rows of titanium staples in 30mm, 45mm, 60mm and 90mm sizes. Cartridges are *usually* sold separately, however we are offering a promo pack for a limited time only. Staple cartridges are available for both standard tissue (3.5mm) and thick (4.8 mm) and Stapler can be cold sterilised and reused.

<u>GVPME-MALS30PK</u>	MIRUS Reloadable Linear Stapler - 30mm INCLUDES 3.5mm and 4.8mm Staple Cartridges	Normally \$522 NOW \$400
<u>GVPME-MALS45PK</u>	MIRUS Reloadable Linear Stapler - 45mm INCLUDES 3.5mm and 4.8mm Staple Cartridges	Normally \$522 NOW \$400
<u>GVPME-MALS60PK</u>	MIRUS Reloadable Linear Stapler - 60mm INCLUDES 3.5mm and 4.8mm Staple Cartridges	Normally \$522 NOW \$400
<u>GVPME-MALS90PK</u>	MIRUS Reloadable Linear Stapler - 90mm INCLUDES 3.5mm and 4.8mm Staple Cartridges	Normally \$522 NOW \$400

(Available until stocks last | Excludes GST and Freight)

GVP
General Vet Products

EOFY SALE!



MIRUSTM AMAZING Linear Cutter Packs Package Deal

MIRUSTM - The Symbol of Excellence which ensures superior surgical outcomes.

- MIRUS Linear Cutter and Reloads are designed for cutting, transecting and stapling of the alimentary canal, linings, tissues and organs. It offers two double staggered rows of titanium staples and simultaneously cuts and staples. MIRUS Linear Cutter & Reloads are available in 60mm, 80mm & 100mm length with two different sizes, 3.8mm and 4.8mm to accommodate various tissue thickness.

<u>GVPME-MLC60PK</u>	MIRUS Linear Cutter - 60mm PLUS <u>Linear CARTRIDGE - 3.8mm and 4.8mm</u>	Normally \$589 NOW \$450
<u>GVPME-MLC80PK</u>	MIRUS Linear Cutter - 80mm PLUS <u>Linear CARTRIDGE - 3.8mm and 4.8mm</u>	Normally \$589 NOW \$450
<u>GVPME-MLC100</u>	MIRUS Linear Cutter - 100mm PLUS <u>Linear CARTRIDGE - 3.8mm and 4.8mm</u>	Normally \$589 NOW \$450

(Available until stocks last | Excludes GST and Freight)

GVP
General Vet Products

EOFY SALE!!

Time to restock your shelves



20% OFF FLEXICARE Products

(Available until stocks last | Excludes GST and Freight)

Circuits . ET Tubes . Adaptors . Rebreathing Bags

GVP
General Vet Products

EOFY

SALE!!

15%
OFF

GVP Stainless Steel Cage orders

Top grade Stainless Steel – SUS304 providing long service life.

Radius turned corners and edges (rounded) eliminating difficult to clean edges and corners.

Single hand locking system making it easy to handle a dog and unlock/lock the door at the same time.

Crimped flooring to help eliminate 'wobble board' sounds.

Easily removable divides.

Featuring our industry leading **Noise dampening system** – see video - You Tube



(Available until stocks last | Excludes GST and Freight)


GVP
General Vet Products



200% OFF



GVP Stainless Steel Tables



(Available until stocks last | Excludes GST and Freight)

COLOURFUL Ecohesive

Bandages

Order a carton
& SAVE

50%



Ecohesive is a breathable, flexible bandage ideal for use by all healthcare professionals and patients with sensitivity to natural rubber latex.

- Available in a wide variety of colours.
- Stretches to approximately 4.5 m in length.
- Self-adherent, cohesive wrap securely sticks to itself and not to skin or other materials.
- Easily tears by hand for custom lengths. No need for scissors.
- Porous construction allows ventilation of the wrapped area.
- Retains dressings without tension around digits or limbs.
- Highly conformable and is simple and painless to remove.
- Non-sterile - Each bandage is individually wrapped.
- Good for 1x autoclave session to use for orthopaedics.



1. GVPTE-EB-FAVE-MP/b - WAS \$556.00
2. GVPTE-EB04-MC-BK/b - WAS \$354.60
3. GVPTE-EB02-MC-BK/b - WAS \$618.00
4. GVPTE-EB02-P-BK/b - WAS \$645.60
5. GVPTE-EB03-MC-BK/b - WAS \$428.80
6. GVPTE-EB03-P-BK/b - WAS \$466.40
7. GVPTE-EB-MLA-MIX/b - WAS \$517.37
8. GVPTE-EB-SA-MIX/b - WAS \$517.37
9. GVPTE-EB-FAVE-MC/b - WAS \$523.40

- NOW \$278.00**
NOW \$177.30
NOW \$309.00
NOW \$322.80
NOW \$214.40
NOW \$233.20
NOW \$258.69
NOW \$258.69
NOW \$261.70