

GVP

General Vet Products



Best Sellers

2019 - 2020

1300 VET GET


www.generalvetproducts.com.au

Log in today to receive
 better-than-catalogue prices
 when shopping online!
www.generalvetproducts.com.au



Contents

Microchips & Scanners	2
Parasite Control	4
IV Therapy	5
Electrosurgery	6
Dispensing	7
Monitoring	8
Tonometry	9
Otosopes and Ophthalmoscopes	10
Bandaging	12
Cages	13
Operating Tables	14
Lights	15



Microchips and Scanners

Home / Microchips & Scanners



Cybortra Microchips

PACK OF 20

The Cybortra Minichip is a one-piece unit that fits easily into the hand and is ergonomically designed for comfortable and efficient use.

The red sleeve must be removed before the applicator stylet is depressed. This avoids inadvertent microchip ejection and also makes identification of used microchip applicators easy, especially when microchipping a litter of pups.

- ISO FDX-B type microchip (required by law for implanting in Australia)
- Manufactured to specifications ISO 11784 and ISO 11785
- Scans at the frequency of 134.2khz
- Each chip has a unique 15-digit microchip number
- Approved for Australia-wide implanting
- Microchip is encapsulated in an inert bioglass
- Diameter 1.4mm, length 8mm
- Individual sterile packaging
- 6 duplicate barcodes included for each microchip

IMPORTANT REGISTRATION INFORMATION

For implanters outside NSW, by law, microchip registration (paper or online) needs to be purchased at the same time, with either Central Animal Records (CAR) or Australasian Animal Register (AAR) – the most used and well-known animal registries in Australia.



New South Wales

Microchips only — no registration required

CODE	PRICE
GVP-CYBORTRA MINI1	\$ 90.00
GVP-CYBORTRA STANDARD1	\$ 82.50

All other states and territories

Microchips require online or paper registration with CAR or AAR

CODE	PRICE
GVP-CYBORTRA AARREGO1	\$ 240.00
GVP-CYBORTRA AARPAPER1	\$ 260.00
GVP-CYBORTRA STANDARDAARREGO1	\$ 220.08
GVP-CYBORTRA STANDARDAARPAPER1	\$ 238.42
GVP-CYBORTRA-CARREGO1	\$ 210.00
GVP-CYBORTRA-CARPAPER1	\$ 260.00
GVP-CYBORTRA STANDARDCARREGO1	\$ 192.57
GVP-CYBORTRA STANDARDCARPAPER1	\$ 238.42

THE MORE YOU BUY THE MORE YOU SAVE

LOG INTO YOUR ONLINE ACCOUNT TO SEE HOW MUCH YOU CAN SAVE!



Portable ISO Microchip Scanner

CODE	PRICE	BULK (3)
GVPDW-PT160	\$ 119.95	\$ 107.96

Used for scanning ISO Microchips (RFID tags) that use the 134.2kHz frequency. It has an OLED display and rechargeable battery. The PT160 can read many kinds of RFID tags, with stable, safe and reliable performance.

Specifications:

- Size: 16 x 9 x 4 cm
- Frequency: 134.2kHz
- Reading distance: Up to 3cm
- Reading standard: ISO 11784/5 FDX-B/ID64
- Charging method: USB
- Charging voltage: 5V
- Charging time: 4-5 hours
- Operation temperature: -15°C to 45°C
- Certification: CE ROHS

Parasite Control

🏠 / Parasite Control / Worming



Fenpral Allwormer Dog 10kg

100 TABLET BOTTLE

CODE	PRICE	BULK (6)
FENTD	\$ 98.50	\$ 86.00

Take 1 tablet every 3 months for the control of all gastrointestinal worms in dogs including:

- Roundworm
- Hookworm
- Whipworm
- Tapeworm
- Hydatid Tapeworm



Fenpral Allwormer Dog 40kg

50 TABLET BOTTLE

CODE	PRICE	BULK (6)
FENTD40	\$ 198.50	\$ 175.00

Take 1 tablet every 3 months for the control of all gastrointestinal worms in dogs including:

- Roundworm
- Hookworm
- Whipworm
- Tapeworm
- Hydatid Tapeworm



Fenpral Allwormer Cat 5kg

100 TABLET BOTTLE

CODE	PRICE	BULK (6)
FENTC	\$ 168.50	\$ 150.00

Take 1 tablet every 3 months for the control of all gastrointestinal worms in cats including:

- Roundworm
- Hookworm
- Tapeworm

Fenpral Dispensing Envelopes

FOR DOG OR CAT, PACK OF 20

CODE	PRICE
FENENV-DOG	\$ 1.80
FENENV-CAT	\$ 1.80

For distributing small quantities of worming tablets.



IV Therapy

🏠 / IV Therapy / IV Pumps



800I

IV INFUSION PUMP

CODE	PRICE
GVPPM-800I	\$ 825.00

The 800I IV Infusion Pump is the latest infusion pump from General Vet Products that brings an easy to use yet sophisticated infusion pump to the veterinary market. It gives vets and nurses the freedom to use an extensive range of standard giving sets, not locked in to using a particular giving set. The unit is easy to calibrate according to the parameters of any giving set tubing that is being used.

A compact, programmable infusion pump suitable for delivery of IV fluid therapy.

Benefits:

- Small, lightweight and compact
- Patented "Double Pulse" compensation technique ensures good linearity even at low delivery rates
- Double CPU system makes the running status super safe.
- Delivery will be cut off, even if control circuit is powered off, and machine will alarm
- Rechargeable back-up battery supports at least 4 hours continuous operation at 25mL/hr
- Different fluid delivery modes available:
 - Fluid delivery rate
 - Delivery time rate
 - Dripping rate
 - 0.1-1500mL/hr flow rate range
- Historical record function, easy to recall operation data
- Alarm for door open, low battery, occlusion, end program, malfunction, etc
- Easy to use rear clamp to attach to IV pole mounts
- Easy to read display with large number format — great for night-time

Dimensions:

174mm(W) x 126mm(D) x 215mm(H)

Weight:

Approx 2.5kg

The 800I IV Infusion Pump works with 10-30 drops/mL giving sets and with varying hardness of tubes from soft to hard. This is a very easy IV Pump to use and calibrate yourself, giving you the flexibility to use whatever giving sets you want to.

Send your giving sets to General Vet Products and we will calibrate the 800I accordingly. First giving set calibration is free. Additional giving sets are calibrated at \$50.00 + GST each.

Note – pricing may vary to the published price from time to time and without notice.

Electrosurgery

🏠 / Equipment / Electrosurgery & Consumables



FT10

ELECTROSURGICAL UNIT

CODE	PRICE
GVPM-VALEC-FT10	\$ 15,985.00
GVPM-VALEC-FT10 EX-DEMO	\$ 8,345.00

The ValleyLab FT10 Electrosurgical Unit (or Energy Platform) is the next generation advanced energy platform with improved LigaSure™ vessel sealing technology and expanded electrosurgical features. Its use spans across all specialties that utilize energy-based devices — inclusive of GYN, Colorectal, Bariatric, General, Urologic and ENT.



Rocker Switch Pencil

REUSABLE — 50 USES

CODE	PRICE
GVPM-E2100	\$ 159.95

Rocker switch pencil with reusable stainless steel hex-locking blade electrode. The elastomeric surface improves the surgeon's grip even when the pencil surface is wet.

- Cord Length: 15' (4.6m) silicone cord
- Accepts all standard 3/32" (2.4mm) shaft electrodes
- Non-sterile
- Steam or ETO sterilizable

Dispensing

🏠 / Medications & Dispensing Containers / Dispensing Containers



Safeseal Medication Bottles

LEAKPROOF & CHILDPROOF

CODE	PRICE	BULK (5)
SS30ML/B	\$ 38.90	\$ 36.96
SS100ML/B	\$ 39.90	\$ 37.91

Safeseal allows you to dispense liquid medications in a leakproof and childproof bottle. The seal is pushed onto the opening of the bottle once you have filled it with your medication to dispense.

Now your client can withdraw accurate and non messy doses using just a syringe inserted into the self closing seal in the top of the Safeseal bottle.

Sold in boxes of 20 including:

- 20 x 30mL or 100mL amber heavy-duty plastic bottles
- 20 x Safeseals
- 20 x childproof caps

Monitoring

Home / Anaesthesia / Monitors / Multi Parameter Monitors



Meditech M3T, M3C & M3E

MINI VET MONITORS

CODE	PRICE	SpO ₂	NIBP	TEMP	EtCO ₂
GVPMT-M3T	\$ 695.00	✓	✓	✓	✗
GVPMT-M3C-S (SIDESTREAM)	\$ 1,495.00	✓	✗	✗	✓
GVPMT-M3C-M (MAINSTREAM)	\$ 1,648.00	✓	✗	✗	✓
GVPMT-M3E-S (SIDESTREAM)	\$ 2,148.00	✓	✓	✗	✓
GVPMT-M3E-M (MAINSTREAM)	\$ 2,148.00	✓	✓	✗	✓

One of the smallest and best presented handheld vital sign monitors on the market.

Features:

- 3.5" high resolution color TFT screen
- NIBP measurement accuracy: ± mmHg
- Pulse-rate range: 40-240bpm
- SpO₂ saturation accuracy range: 70-100%
- Temperature measurement range: 25-45°C, Accuracy: 0.1°C
- EtCO₂ measurement range: 0-15mmHg (0-19.7%) (0-20kPa)
- Power Type: Lithium ion battery, up to 8 hours for normal operation

Includes:

- Reusable SpO₂ sensor with 2 clips
- 1 set of disposable cuffs
- NIBP extension cable
- Reusable temperature probe (GVPMT-M3T only)
- Charging stand charging adapter Uses rechargeable Li-ion battery or AA batteries
- USB connection cable
- User manual
- EtCO₂ adapter (EtCO₂ models only)



Meditech M3S Mini Vet

CODE	PRICE
GVPMT-M3S	\$ 448 .00

Features:

- 3.5" high resolution colour TFT screen
- Reusable SpO₂ sensor with 2 clips
- Charging stand and charging adapter
- Uses rechargeable Li-ion/AA batteries for 8 hours of monitoring



Tonometry

Home / Equipment / Tonometers



TONOVET

LEAKPROOF & CHILDPROOF

QUANTITY	CODE	PRICE
	GVPS-TONOVETPLUS	\$ 4,145.00
Box of 100	GVPS-TONOVETPROBES	\$ 171.85



The Icare TonoVet utilises a technology called 'rebound' tonometry. The portable hand held unit is about 1/4" away from the cornea. At the push of a button, a tiny probe painlessly and rapidly reaches out and touches the cornea. Six rapid measurements are taken within seconds. Most of the time the animal will not even react and certainly no local anaesthesia is required. An audible beep tells of a successful result, which is displayed digitally.

The Icare TonoVet has been used on all veterinary patients including equine. The TonoVet is self-calibrating with its own software.

Features:

- Easy to use and fully automated
- Provides faster results
- More comfortable for the patient, no topical anaesthetic required
- Higher accuracy at higher and lower IOP
- Instructional DVD included
- Held upright in a natural, ergonomic position like an ophthalmoscope versus horizontally like a flashlight



Specifications:

- Weight: 155g without batteries
- 250g with batteries
- Power Supply: 4 x AA batteries (approx. 1,000 measurements per battery set)
- Measurement Range: 1-99 mmHg
- Size: 13-32mm (W) x 45-80mm (D) x 230mm (L)

Kit includes:

- Tonometer
- Padded storage case
- Batteries
- Probes (100pk)
- 2-year limited warranty



Otoscopies & Ophthalmoscopes

Home / Equipment / Otoscopes & Ophthalmoscopes / Welch Allyn & Compatible



Diagnostic Set with USB Charger

CODE	PRICE
GVPAZ-468E4200	\$ 736.25

Includes:

- USB charging handle
- LED coaxial ophthalmoscope
- LED operating otoscope
- Hard case



Diagnostic Set with Single Pod

CODE	PRICE
GVPAZ-468E5201/b	\$ 796.00

Includes:

- Single charging pod
- LED operating otoscope
- Lithium rechargeable handle
- LED coaxial ophthalmoscope
- Hard case

Operating Oscope Head

LED 3.5V



CODE	PRICE
GVPAZ-462E0010	\$ 290.50

Compatible with GVPZ-4890010 and GVPZ-4900010, and other leading brands of 3.5V otoscope handles.

Ophthalmoscope Pro Coaxial

LED 3.5V



CODE	PRICE
GVPAZ-411E0010	\$ 343.00

Compatible with Zumax and other 3.5V ophthalmoscope handles.

Otoscope Handle

WITH USB CHARGER



CODE	PRICE
GVPAZ-4900015	\$ 345.80

- Metal handle
- Rechargeable lithium ion battery

Otoscope Handle

FOR DESK POD



CODE	PRICE
GVPAZ-4900010	\$ 175.68

- Metal handle
- Rechargeable lithium ion battery

Lightweight Oscope Handle

RECHARGEABLE 3.5V



CODE	PRICE
GVPAZ-4890010	\$ 118.00

- Desktop pod charging only
- Compatible with Welch Allyn 3.5V heads

Wall-mount Handle Set



CODE	PRICE
GVPAZ-485-0100	\$ 396.00

- Double wall-mount otoscope/ ophthalmoscope handle set - 3.5V
- Compatible with Zumax and other 3.5V heads

AS USED IN NSW PUBLIC HOSPITAL A&E

COMPATIBLE WITH RECENT WELCH ALLYN™ MODELS



Desk Chargers

FOR ZUMAX OTOSCOPES AND OPHTHALMOSCOPES

	CODE	PRICE
SINGLE	GVPAZ-4871000	\$ 150.15
DOUBLE	GVPAZ-4881000	\$ 364.60

- 12V
- For use with handles GVPZ-4900010 & GVPZ-4890010, also compatible with handles from other leading brands

Power Cable

FOR DESK CHARGERS



CODE	PRICE
GVPAZ-9801229	\$ 77.06

- 12V
- For chargers GVPZ-4871000 and GVPZ-4881000

Lithium Ion Battery

FOR 3.5V LITHIUM HANDLES



CODE	PRICE
GVPAZ-4900980	\$ 84.50

Ophthalmoscope Globe

P4900, 3.5V



	CODE	PRICE
Halogen	GVPAZ-4100910	\$ 36.26
LED	GVPAZ-410E0910	\$ 105.00

Otoscope Globe

T3100, 3.5V



	CODE	PRICE
Halogen	GVPAZ-4610910	\$ 29.00
LED	GVPAZ-461E0910	\$ 78.50

Accessories

LED, 3.5V



	CODE	PRICE
Finoff Illuminator	GVPAZ-472E0010	\$ 156.80
Nasal Illuminator	GVPAZ-473E0010	\$ 148.50
Nasal Bivalve Speculum	GVPAZ-474R0010	\$ 187.50

Compatible with Zumax and other 3.5V handles.

Operating Tables

Home / Equipment / Tables - Operating



Vertical Lifting Operating Table

WITH V-TOP MECHANISM

CODE	PRICE
GVPTE-FT-886	\$ 3,495.00

V-top Operation Table - powered by a sturdy and reliable Danish lift column.

Features:

- A two-panel table-top system that adjusts from horizontal to near vertical to hold the pets in the desired position during surgery
- Overall size of the top: 140cm x 50cm
- Powerful and durable electric lifting actuator holds 120kg (240lbs) on the surgery top
- Two sliding, quick release easy controls are provided on each side of the rail to offer secure pet restraint
- 10cm (4") stainless steel trough encourages fluid run off into a drain
- Table can be lifted from 61cm to 101cm
- Worktop Swivel Angle: 0-45°

All prices exclude GST, freight charges may apply, see website for Terms & Conditions

Lights

Home / Lights / Operating



Burton Aim LED Surgery Light

burton

AVAILABLE IN 4 CONFIGURATIONS

CONFIGURATION	CODE	PRICE
Single Ceiling	GVPB-AIMLEDSC	\$ 5,692.40
Dual Ceiling	GVPB-AIMLEDDC	\$ 11,289.00
Wall Mount	GVPB-AIMLEDWM	\$ 5,692.40
Floor Stand	GVPB-AIMLEDFL	\$ 5,861.50


The AIM LED uses the latest advances in state-of-the-art lighting technology to achieve incredible performance and unsurpassed value and energy efficiency. Using 70% less energy than traditional light source products and with 10 times greater lamp life (20,000 hours), you will save more than \$1,000 in bulb replacements alone.

- Illuminance: 45,000lux (4181fc) at 1m
- Color temperature: 4300K
CRI (Color Rendering Index): 92
- Diameter of lighthouse: 20" (51cm)
- Light field diameter: 13" - 15" (33cm-38cm)
- Depth of illumination: 54" (137cm)
- Focusing: Adjustable using center handle
- Number of LED modules: 3
- Light sources: Light Emitting Diodes
- (LED) Rated life of LED lamp: 20,000 hours
- Power: 84 Watts
- Yes Certifications/approvals: UL60601-1 / IEC 60601-1 / IEC 60601-1-2 / IEC 60601-2-41 CAN/CSA C22.2 601-1 M90 CE marked
- Product warranty: 5-year parts and labor
Origin: Assembled in USA

All prices exclude GST, freight charges may apply, see website for Terms & Conditions

CONTACT US

Mail	PO Box 417 Fairy Meadow 2519 NSW
Phone	1300 838 438 / 1300-VET-GET
Phone	02 4229 1700
Email	sales@generalvetproducts.com.au
Web	www.generalvetproducts.com.au
LinkedIn	au.linkedin.com/company/gvp-pty-ltd



GVP
General Vet Products

Cover image: Helena Lopes
@wildlittlethingsphoto