

GVP

General Vet Products



IV Therapy

PUMPS AND ACCESSORIES

2019 - 2020



1300 VET GET

www.generalvetproducts.com.au



Contents

IV Pumps	3
Administration Sets	8
IV Pump Accessories	12
IV Cannulas	13
Ports, Plugs & Stoppers	14

**LOG IN TODAY TO RECEIVE
BETTER-THAN-CATALOGUE PRICES
WHEN SHOPPING ONLINE!
www.generalvetproducts.com.au**

SPRING SALE!



CHECK OUT OUR
BRAND NEW
INFUSION PUMPS

GET ONE TODAY FOR JUST \$825*

BUY NOW AND RECEIVE:

**20 INTRAFIX PRIMELINE GIVING SETS
+ 2 BOXES OF UNIVERSAL EXTENSION SETS**

**OVER \$100 OF VALUE! **
ABSOLUTELY FREE!!**

**PLUS!
WE WILL
CALIBRATE
YOUR PUMP**
TO SUIT THE GIVING SETS
AT NO EXTRA CHARGE!

OFFER VALID UNTIL 31ST DECEMBER 2019

*All prices exclude GST ~ Freight charges may apply, please see website for Terms & Conditions
**Only 1 package consisting of: 20 x giving sets & 2 x boxes extension sets (80pieces) per practice

Infusion Pumps



8001
IV INFUSION PUMP

CODE	PRICE
GVPPM-8001	\$ 825.00

The 8001 IV Infusion Pump is the latest infusion pump from General Vet Products that brings an easy to use yet sophisticated infusion pump to the veterinary market. It gives vets and nurses the freedom to use an extensive range of standard Giving Sets, not locked in to using a particular giving set. The unit is easy to calibrate according to the parameters of any Giving Set tubing that is being used.

A compact, programmable infusion pump suitable for delivery of IV fluid therapy.

Benefits:

- Small, lightweight and compact
- Patented "Double Pulse" compensation technique ensures good linearity even at low delivery rates
- Double CPU system makes the running status super safe. Delivery will be cut off, even if control circuit is powered off, and machine will alarm
- Rechargeable Back-up Battery supports at least 4 hours continuous operation at 25mL/hr
- Different fluid delivery modes available:
 - Fluid delivery rate
 - Delivery time rate
 - Dripping rate
- 0.1-1500mL/hr flow rate range
- Historical record function, easy to recall operation data
- Alarm for door open, low battery, occlusion, end program, malfunction, etc
- Easy to use rear clamp to attach to IV pole mounts
- Easy to read display with large number format - great for night-time

Dimensions:
174mm(W) x 126mm(D) x 215mm(H)
Weight:
Approx 2.5kg

The 8001 IV Infusion Pump works with 10-30 drops/mL Giving Sets and with varying hardness of tubes from Soft to Hard. This is a very easy IV Pump to use and calibrate yourself, giving you the flexibility to use whatever Giving Sets you want to.

Send your giving sets to General Vet Products and we will calibrate the 8001 accordingly. First Giving Set calibration is free. Additional Giving Sets are calibrated at \$50+GST each.

Note – pricing may vary to the published price from time to time and without notice.

* All prices exclude GST ~ Freight charges may apply, please see website for Terms & Conditions

Infusion Pumps



Niki V4

IV INFUSION PUMP

CODE	PRICE
GVPE-5570(a)	\$ 1,595.00

A compact, programmable infusion pump suitable for delivery of fluids, including total parenteral nutrition (TPN) and intravenous fluid therapy. Niki V4 Volumetric Infusion Pump uses standard 10 drops/mL IV Giving Sets.

Benefits:

- Small, lightweight and compact
- Macro and Micro infusion capability
- Rechargeable internal battery
- Up to 10 programme settings
- Niki V4 Volumetric Infusion Pump uses standard 10 drops/mL giving sets

Infusion programs:

- Micro
- Macro
- 10-Steps
- TPN

Money-Saving Rechargeable Battery

A recent cost analysis estimates a savings of \$600.00 per pump per year using the rechargeable battery.

Piston Mechanism

Delivers $\pm 5\%$ accuracy and 1.8 μ L linear infusion even at low flow rates using dedicated infusion sets.

Post-Occlusion Bolus Reduction System

Prevents delivery of false boluses following a downstream occlusion event.

Precise In-Line Pressure Detection

Detects downstream and upstream pressure.

Note – pricing may vary to the published price from time to time and without notice.

Dimensions:

138mm(W) x 138mm (D) x 72mm (H)

Weight:

Approx 1.1kg

Do you have an online account with us?

We offer great online discounts and promotions to our customers when shopping on our website –

www.generalvetproducts.com.au

Log in or Sign Up today!

Syringe Pumps



SK-5000

SYRINGE PUMP

CODE	PRICE
GVPBL-DE10-2500	\$ 790.00

Built-in double CPU

Each CPU can serve as backup for the other, guaranteeing the stability of the system and the safety of patients even when one CPU breaks down.

Powerful alarm system with both audible and visible alarm signs

- Embedded pressure sensor for occlusion detection
- Multiple alarm signs to alert operators

Extreme flexibility in the use of syringes

Almost all brands of standard syringes can be applied and automatically recognized by the module to save your cost and time.

Wireless linkage with infusion supervision system

Greatly reduces workload and turnaround time, and ensures the safety of patients.

Seven-item protection design

Water-proof, moisture-proof, anti-interference, anti-shock, anti-collision, anti-corrosion, anti-high and low temperature.

Wide-input power supply

Suitable for a broad range of circumstances and countries.

Compatible with: 5, 10, 20, 30, 50 & 60mL syringes.

VTBI: from 0.1 to 999.9mL

Accuracy: +/-2%

Bolus Rate: 100ml/h to 500mL/h

KVO rate: 0.5 ml/h

Automatic detection of syringe type.

Dimensions:

310mm (W) x 125mm (D) x 115mm (H)

Weight:

< 3kg

Occlusion:

Low: 300 +/- 125 mmHg

Medium: 525 +/- 125 mmHg

High: 900 +/- 180 mmHg

Flow rate range:

5mL syringe: 0.1 - 150 mL/h

10mL syringe: 0.1 - 300 mL/h

20mL syringe: 0.1 - 600 mL/h

30mL syringe: 0.1 - 900 mL/h

50 and 60mL syringe: 0.1 - 1,500 mL/h

Display

Screen: 2.0 inch B/W LCD

Information: Rate, accumulated injection volume, syringe size, battery capacity, bed no., AC power indicator, battery charge indicator

Battery

Type: 7.4V, 1600mAh lithium polymer battery

Life: More than 2 hours' continuous running (5 ml/h); rechargeable

Syringe Pumps

Anitek M200A SYRINGE PUMP



CODE	PRICE
GVPE-5494	\$ 1,630.00

The Anitek M200A Syringe Pump is a fixed capacity pump, which delivers drugs in a specific time, without being effected by resistance inside the infusion tube.

Equipped with a micro-propulsion system, it can achieve a wide speed regulation range, through the control circuit.

When any syringe size (5mL, 10mL, 20mL, 30mL and 50mL) is loaded onto the pump, the identification system can automatically identify the specification of the loaded syringe.

At this moment, the user only has to set the required speed and press the ON/OFF button to start the pump.

The Anitek M200A Syringe Pump has the following features:

- Suitable for 5mL, 10mL, 20mL, 30mL and 50mL syringes
- Precise pressure monitoring system displays real time, dynamic changes in the syringe line pressure
- Reliable pushrod combination, ensures safe operation of the syringe and avoids potential dangers
- Clamp, on rear of unit for attaching to pole
- Portable design makes it convenient and sturdy

Dimensions:

308mm (W) x 144mm (W) x 156mm (H)

Weight:

Approx 2.2kg

Speed:

- 50ml syringe: 0.1~1200 mL/hr (0.1 mL/h each level)
- 30ml syringe: 0.1~600 mL/hr (0.1 mL/h each level)
- 20ml syringe: 0.1~400 mL/hr (0.1 mL/h each level)
- 10ml syringe: 0.1~300 mL/hr (0.1 mL/h each level)
- 5ml syringe: 0.1~100 mL/hr (0.1 mL/h each level)

Accuracy: ±3%

Fast infusion:

- 1200 mL/hr (50mL syringe)
- 600 mL/hr (30mL syringe)
- 400 mL/hr (20mL syringe)
- 300 mL/hr (10mL syringe)
- 100 mL/hr (5mL syringe)

Total amount:

0.1mL~9999.9mL
(0.1mL each level)

Alarm:

Residual alarm, completed injection alarm, tube blockage alarm, unclamped cylinder alarm, improperly installed putter alarm, out-of-limit speed alarm, alarm of infusion equivalent to limited amount, disconnected power line alarm, battery power depletion alarm, system error alarm, neglected operation alarm.

Power supply:

AC220V, 50Hz; in-built DC 11.1V rechargeable lithium battery; the fully charged battery can sustain the pump operation for more than 3 hours at the speed of 5mL/h (specified by EN 60601-224:1998 standard).

*All prices exclude GST ~ Freight charges may apply, please see website for Terms & Conditions

Fluid Warmers

i-Warm™ IV FLUID WARMER



CODE	PRICE
GVPE-5535-5	\$ 395.00

The iWarm™ IV Fluid Warmer is an economical fluid warmer with temperature display for standard IV sets.

Features:

- Portable and lightweight
- Convenient hanging belt to hang on IV pole
- Current temperature displayed in degrees Fahrenheit or Celsius
- Temperature overheating protection: indicators for warming (orange), power (green) and overheating (red)
- Reduces prep time required to warm fluid
- Can monitor temperature second to second
- Reduces patients' discomfort due to cold fluid
- Compatible with standard IV sets
- Quickly brings the temperature level of the fluid up to that of the patient
- Not recommended for use with blood products, and not suitable for flow rates less than 50mL/h
- Backup overheating protection with power interrupted at 45°C with visual and audio alarm.
- 1-year warranty



Specifications:

- IV tubing size: 3.0 - 5.0mm
- Flow rate range: 50-500mL/h
- Power requirement: 240V AC 60Hz
- Weight w/power cord: 0.65kg
- i-Warmer dimension: 19 x 9 x 6cm
- Suggested operating temp: 20-25°C
- Heating method: Dry heat exchange

If you cannot find the products that you need, then please don't hesitate to contact our Sales team today.

General Vet Products has always been proud of being able to source almost any product for you.

*All prices exclude GST ~ Freight charges may apply, please see website for Terms & Conditions

Giving Sets



B | BRAUN Intrafix® Primeline GIVING SET 180cm

THAT'S JUST \$1.33 PER SET!

CODE	PRICE per 10
GVPBA-4062181	\$ 13.27

Product Features:

- Large 15µm filter, provides high flow rate and protects against larger particles
- Advanced precise roller clamp - good to control flow rate, comfortable to hold. Spike dock for secure spike discard. Housing colour "white"
- Screw tight Luer Lock fitting - prevents accidental disconnection
- Soft, clear PVC tubing with standard lumen 3mm (ID) x 4.1mm (OD)
- Standard length 180cm, DEHP free plasticiser with superior, well documented safety properties
- Ergonomically designed 3-part drip chamber
- Universal, vented piercing spike – spike IV bags with ease and prevent contamination
- Bacteria-tight proven air vent, closable with handy snap cap
- Fully transparent upper allows quick drop rate check
- Push ring for comfortable and secure gripping – ideal for attachment of Infusomat® pump drip sensor
- Extended, detached drip element provides even drop formation; 20 drops = 1mL ± 0.1mL (aqua dist)
- Elastic pump chamber enables one-click fluid level setting
- Offers a high fluid level as guard against air entrainment
- Needle based, air-free latex Eurofix® injection site with handy safe grip plate
- All design features in accordance with International

BBraun Intrafix® Primeline Giving Sets are a Standard 20 drops/mL IV Set for administering all kinds of infusions. They are a universal IV administration set with air vent for gravity infusions and pressure infusions with non-dedicated infusion pumps (pressure resistant up to 2 bar).

Needle based, air-free latex Eurofix® injection site with handy safe grip plate.

* All prices exclude GST ~ Freight charges may apply, please see website for Terms & Conditions

Giving Sets



B | BRAUN Sangofix ES® GIVING SET 180cm

THAT'S JUST \$2.32 PER SET!

CODE	PRICE per 10
GVPBA-4100514	\$ 23.23

Manufactured for the intravenous administration of blood and blood components but can also be used for the administration of other physiological fluids eg saline, Hartmanns etc.

The drip chamber is equipped with a standard blood filter suitable for blood and blood components. Removes blood clots greater than 200µm. Gravity administration set suitable for 3 to 4 blood units. 20 drops/mL.

'Prime Stop' Function allows priming of the fluid line and expulsion of air without draining fluid from the giving set end. Once primed the green 'Prime Stop' is removed and discarded.

With standard lumen 3 x 4.1mm, standard length 180cm, DEHP free plasticiser with superior, well-documented safety properties.

Niki IV Giving Set

COMPATIBLE WITH V4 FLUID PUMPS (50s)

CODE	PRICE per 50
GVPE-5555	\$ 120.00

- Standard Administration Set, suitable for use with Niki V4 Fluid Pumps
- 170cm with injection port and male luer slip
- 10 drops per ml
- Sold in boxes of 50
- Priced at \$2.40 per set



* All prices exclude GST ~ Freight charges may apply, please see website for Terms & Conditions

Extension Sets



Universal

EXTENSION SETS

LESS THAN \$1 PER SET!

CODE	LENGTH	PRICE per 40
GVP-GEXT160U	160cm	\$ 39.50
GVP-GEXT90U	90cm	\$ 39.50

Our universal extension sets are everything you want on an extension set with dual connection types at the catheter joining end (male luer slip and luer lock).

- Common bore extension set
- Comes in 160cm or 90cm length
- Injection ports joined solidly in place – do not come apart unlike some other brands
- Clip lock
- Sold in boxes of 40

Extension Sets



Universal

MINIMUM VOLUME



CODE	LENGTH	PRICE per 40
GVP-GEXT160	160cm	\$ 53.50

Our universal minimum volume extension sets have the following features:

- Ideal for low fluid rates with small animal patients
- Highly flexible, kink resistant tubing
- Dual connection types at catheter joining end - male luer slip and luer lock
- Comes in 160cm length
- Injection ports joined solidly in place - do not come apart unlike some other brands
- Clip lock
- Sold in boxes of 40

THAT'S JUST \$1.34 PER SET!

B | BRAUN Heidelberg

EXTENSION SETS

CODE	LENGTH	PRICE per 100
GVPBA-4033809	22cm	\$ 203.69
GVPBA-4034589	70cm	\$ 416.08
GVPBA-4055128	75cm	\$ 229.78
GVPBA-4055136	140cm	\$ 273.24

BBraun Heidelberg extension tubing is a common IV giving set that can be extended easily without length limitations.

- Excellent quality tubing at an affordable price
- Luer lock connection
- No limited flow-rate for all standard applications
- Available in lengths: 22cm, 70cm, 75cm and 140cm
- Sold in boxes of 100

* Should not be used as extension line for syringe pumps



* All prices exclude GST ~ Freight charges may apply, please see website for Terms & Conditions

B | BRAUN Heidelberg

MINIMUM VOLUME

CODE	LENGTH	PRICE per 100
GVPBA-4097572	180cm	\$ 223.56
GVPBA-9500065	140cm	\$ 240.00
GVPBA-9500057	75cm	\$ 390.00

BBraun minimum volume extension sets only require 1.2mL to prime the line.

- Male luer slip and lock at one end and Female luer slip at the other end.
- Premium product at affordable price in the following lengths: 75cm, 140cm, 180cm
- Sold in boxes of 100 (75cm length sold in boxes of 200)



* All prices exclude GST ~ Freight charges may apply, please see website for Terms & Conditions

IV Pole and Pump Mounts

Freestanding IV Pole

ON 5-STAR BASE, WITH CASTORS

CODE	PRICE
GVPTE-1005	\$ 159.50

Features:

- Sturdy IV pole on 5-star base with castors
- Solid, smooth, high density polyethylene base
- Ramshead shaped hooks able to hold multiple IV bags
- Height adjustable
- Stainless steel poles
- **Available September 2019**

NEW STOCK!



GVPTE-1005

Cage Mount IV Pole

AND PUMP MOUNT

CODE	PRICE
GVPTE-1000	\$ 128.50

Features:

- Sturdy IV pole that hooks and braces to cage door
- Ramshead shaped hooks able to hold multiple IV bags
- Height adjustable
- Stainless steel poles



GVPTE-1000

IV Cannulas

BRAUN Intracan Certo

NON-WINGED CANNULAS



Technical Specifications

PUR (Polyurethane), Certo®

- Ergonomic design for one-handed insertion
- High stability during initial insertion and high flexibility inside the vein
- Sharp backcut bevel for minimal puncture trauma
- Transparent segments provide confirmation of successful catheter placement when withdrawing the needle
- High adjustment to vessel anatomy and thereby highest comfort for the patient
- Fast and easy identification of gauge size even in low lighting conditions
- Latex and PVC free
- Removable rear plug can be inserted to prevent cannula blood seep
- Sold in boxes of 50



SIZE GUIDE

Gauge 14 and 16

Large Animal
Use for fluid administration and IV access

- Major surgical procedure
- Large volume fluid or blood replacement
- Rapid, bolus infusions

Gauge 18

Small Animal – medium to large breed dogs
Use for large volume and rapid fluid administration

- General surgical procedures
- Transfusion of blood and blood products

Gauge 20 and 22

Small Animal – small to medium breed dogs
Use for IV fluid administration

- Appropriate for most infusion therapies

Gauge 24

Small Animal – cats and pocket pets

- Used for small, fragile veins

\$1.20 PER CANNULA

CODE	SIZE	PRICE per 50
GVPBA-4251369	14G x 50mm (2")	\$ 59.95
GVPBA-4251350	16G x 50mm (2")	\$ 59.95
GVPBA-4251342	18G x 32mm (1.25")	\$ 59.95
GVPBA-4251334	18G x 45mm (1.75")	\$ 59.95
GVPBA-4251326	20G x 32mm (1.25")	\$ 59.95
GVPBA-4251318	22G x 25mm (1")	\$ 59.95
GVPBA-4251300	24G x 19mm (.75")	\$ 59.95

* All prices exclude GST ~ Freight charges may apply, please see website for Terms & Conditions

* All prices exclude GST ~ Freight charges may apply, please see website for Terms & Conditions

Injection Ports



B | BRAUN

Injection Port

INT STOPPER WITH LUER LOCK

CODE	PRICE per 100
GVPBA-4238010	\$ 30.38

Product Features:

- Proximal self-sealing injection port
- Distal male luer lock connector with recessed luer taper
- Components do not contain DEHP or Latex
- Great for plugging an IV catheter to maintain IV access when you have taken your patient off IV fluids
- Sold in boxes of 100

JUST 30¢ EACH!

Injection Port

INT STOPPER WITH LUER SLIP

CODE	PRICE per 100
GVPBA-BFI10010LSP	\$ 38.90

Product Features:

- Luer slip injection port
- 2.2cm length
- Latex injection port
- Great for plugging an IV catheter to maintain IV access when you have taken your patient off IV fluids
- Sold in boxes of 100

NEW STOCK!
LESS THAN 39¢ EACH!



Stoppers & Connectors

B | BRAUN

Combi Stopper

CLOSING CONES

CODE	PRICE per 100
GVPBA-4495209	\$ 19.80

White Combi Stopper closing cones from BBraun:

- Luer Lock fitting male and female
- Made of PE (Polyethylene)
- Sterile closing of all open IV-accesses in IV-Sets
- Sterile closing of pre-filled syringes
- Sold in boxes of 100



LESS THAN 20¢ EACH!

MILA
INTERNATIONAL, INC.

3-Way Stopcock

CODE	PRICE (each)
GVPML-2400	\$ 7.50
GVPML-2405	\$ 7.50

MILA Three Way Stopcock

- 2400 – Three-way Stopcock with Slip Luer
- 2405 – Three-way Stopcock with Luer Lock
- Sterile
- Sold individually



B | BRAUN

Transofix

TRANSFER SPIKES

CODE	PRICE per 50
GVPBA-4090500	\$ 57.95

Transfer spike for large containers

- Fits into rubber stoppers 6.3mm and bag set ports
- Allows easy flow even with non-evacuated glass vials
- Large internal diameter
- Provides secure connection of two containers
- Large grip flanges
- Sold in boxes of 50



THAT'S JUST \$1.16 EACH!

* All prices exclude GST ~ Freight charges may apply, please see website for Terms & Conditions

* All prices exclude GST ~ Freight charges may apply, please see website for Terms & Conditions

A veterinarian in blue scrubs is holding a brown and white dog. The dog has its tongue out and is looking towards the camera. The veterinarian's hands are visible, gently holding the dog's neck and chest. The background is a light-colored wall with a blue horizontal stripe.

CONTACT US

Mail	PO Box 417 Fairy Meadow 2519 NSW
Phone	1300 838 438 / 1300-VET-GET
Phone	02 4229 1700
Email	sales@generalvetproducts.com.au
Web	www.generalvetproducts.com.au
LinkedIn	au.linkedin.com/company/gvp-pty-ltd



www.generalvetproducts.com.au
1300-VET-DE

GVP

General Vet Products



2019 - 2020